



# Face shield

Regular fit




Extra space




[support.apple.com/faceshield](https://support.apple.com/faceshield)

If you need additional assistance with your face shield, please email [faceshield@apple.com](mailto:faceshield@apple.com).

## ① Align your preferred slots on the face shield and forehead band.



## ② Thread the strap through the slots and pull to adjust.



## ③ Repeat on other side.

**Disclaimer:** Intended to be worn with personal protective equipment per your institution's standards. This product has not been tested or qualified to prevent or reduce infection, and does not provide particulate filtration. Not intended for use in surgical settings. This product has not been disinfected or cleaned. Non-sterile. Avoid exposure to high heat. Should not be shared with others. This product has been evaluated to be compatible with common cleaning chemicals. See website for more details.

This product has not been FDA cleared or approved. The product has been authorized by FDA under an EUA for use by healthcare providers as personal protective equipment. This product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use under Section 564(b)(1) of the Act, 21 USC 360bbb-3(b)(1) unless the authorization is terminated or revoked sooner.